



The EUROSETS™ NEW CHALLENGE

"NEGATIVE PRESSURE WOUND THERAPY"

WATERLILY™ is his name chosen by Eurosets for the new generation of devices dedicated to the negative pressure wound therapy.

The official launch is planned during the EACTS congress (European Association of Cardio-Thoracic Surgery) that will take place in Geneva from 11th to 15th of September.

WATERLILY™ is a dynamic and exclusive non invasive system, which, by means of a negative pressure on the wound, allows to:

- increased the local blood flow circulation
- reduce the edema
- stimulate the formation of granulation tissue
- stimulate the cell proliferation
- reduce bacterial load

The fields of application are chronic ulcers at the lower limbs, diabetic foot ulcers, pressure ulcers, engraftment of skin grafting, surgical dehiscence and traumatic injuries.



Just think that in Italy ulcers affect 1 million to 3 million people, with an overall incidence between 0.4 and 4% in the elderly population.

Besides, 15% of diabetic patients presents lifetime ulcers of the lower limbs and 12-24% of these may require further amputation.

The negative pressure generated by **WATERLILY™** is controlled by a microprocessor and is transmitted through a dedicated and innovative dressing to the wound.

One Reservoir collects exceeding exudates leaving to the physician the complete system control at any time.

The lightness of the device, combined with the advanced dressing, promotes rapid and effective wounds healing, guaranteeing to the Patients during the therapy, a good quality of life.

Full system **WATERLILY™** detail on www.eurosets.it



INNOVATION AND TECHNOLOGY FOR A BETTER PATIENT OUTCOMES

